



Laboratory Specimen Requirements

TEST	McKesson Order Code	Lab Order Code (Lab Use Only)	Container	Amount	Special Consideration
Profiles:					
Basic Metabolic Panel	BMP	BMP	Light Green (Lithium Heparin)		
Comp Metabolic Panel	CMP	CMP	Light Green (Lithium Heparin)		
Electrolyte Panel(EP)	EP	EP	Light Green (Lithium Heparin)		
Hepatic Panel	HFP	HFP	Light Green (Lithium Heparin)		
Hepatitis Panel	HEPAP	HEPP	SST (Serum Separator)		
Lipid Panel	CARE2	LIPID	Light Green (Lithium Heparin)		
Renal Function Panel	RFP	RFP	Light Green (Lithium Heparin)		

Urinalysis	UA	UMACR	Urine		
Urinalysis with reflex	URINEC	URINEC	Urine		Clean catch specimen (Live 3/2/2015)
Hemo/Coag					
APTT	APTT	PTT	Plastic Blue (Buffered NA Citrate)	full tube 2.7ml	
PLATELET FUNCTION ASSAY (PFA)	PFA	PFA	Glass Blue (Buffered Citrate NA)	full tube 4.5ml	2 tubes/ no 23 gauge needle
CBC WITH DIFF	CBCD	CBC	Lavender (EDTA)		
CBCD without Diff	CBC	CBCD	Lavender (EDTA)		
D-Dimer	DDIM2	DDIM2	Plastic Blue (Buffered NA Citrate)	full tube 2.7ml	
Ehrlichia IgG/IgM	EHRL	EHRL	SST (Serum Separator)		
Eosinophil Count	EOCT	EOCT	Lavender (EDTA)		
Fibrinogen	FBG	FIBR	Plastic Blue (Buffered NA Citrate)	full tube 2.7ml	
Prothrombin Count (PTR)	PTR	PTR	Plastic Blue (Buffered NA Citrate)	full tube 2.7ml	
Sed Rate (ESR)	SR	SR	Black		
Blood Parasite Screen	BLDPAR	BLDPAR	Lavender (EDTA)		

TEST	McKesson	Lab	Container	Amount	Special Consideration
LDH	LDH	LDH	Light Green (Lithium Heparin)		
Lithium	LI	LI	Gold (Serum Separator)		
Magnesium	MG	MG	Light Green (Lithium Heparin)		
Microalbumin (MA)	MA/R	MA	Urine		
Phosphrous	PO4	PHOS	Light Green (Lithium Heparin)		
Protein, total	TP	TP	Light Green (Lithium Heparin)		
Pregnancy Test/Qualitative	HCG	HCG	Gold (Serum Separator)		
Pregnancy Urine/Qualitative	HCGU	HCGU	Urine Cup or Yellow Urine Tube		
Pregnancy Test/Quantitative	B-HCG	BHCG	Light Green (Lithium Heparin)		
PSA	PSA	PSA	Light Green (Lithium Heparin)		
PSA Free and Total	PSATF	PSAFT	SST (Serum Separator)		
T3, Free	FT-3	FT3	Gold (Serum Separator)		
T3 Uptake	T3UP	T3U	Light Green (Lithium Heparin)		
T3 Total	T3TOT	F3TOT	Gold (Serum Separator)		
T4	T4	T4	Light Green (Lithium Heparin)		
T4,FREE	FT-4	FT4	Light Green (Lithium Heparin)		
Testosterone Free and Total	TESTF	TESTF	Glass Red (No additive)		
Testosterone Total	TESTO	TESTO	Gold (Serum Separator)		
Theophylline	THE	THEO	Light Green (Lithium Heparin)		
Triglycerides	TG	TRIG	Light Green (Lithium Heparin)		
TSH-3rd generation	TSH2	TSH	Light Green (Lithium Heparin)		
Uric Acid	URIC	URIC	Light Green (Lithium Heparin)		
Vitamin B12	VB12	VB12	Gold (Serum Separator)		
Vitamin B1	VITB1	VITB1	Lavender (EDTA) <u>Protect from Light</u>		Thiamine/Protect from Light
Vitamin D	VITDT	VITDT	Gold (Serum Separator)		
Vitamin D fractions	VIT-D	VITD	SST (Serum Separator)		
Vitamin D 1,25	MISCE	MISCE	SST (Serum Separator)		
Vitamin E	VITE	VE	SST (Serum Separator)		
Immunology					
ANA Screen	ANA	ANAQ	SST (Serum Separator)		
Anti DS DNA	DSDNA	DSDNA	SST (Serum Separator)		
Anti JO-1	AUTJO	AUTJO	SST (Serum Separator)		
Anti SCL-70	SCL70	SCL70	SST (Serum Separator)		

TEST	McKesson	Lab	Container	Amount	Special Consideration
Anti SM/RNP	ANTSM	ASM2	SST (Serum Separator)		
Anti SSA	AUTSA	AUTSA	SST (Serum Separator)		
Anti SSB	AUTSS	AUTSS	SST (Serum Separator)		
ASO Screen	ASO	ASOTT	SST (Serum Separator)		
H Pylori IGG	HYPLO	HYPLO	SST (Serum Separator)		
Lyme with Reflex to WB	LYME	LYME	SST (Serum Separator)		
Lyme Western Blot	LYMES	LYMES	SST (Serum Separator)		
Mono Screen	MONO	MONO	Gold (serum Separator)		
Rheumatoid Factor	RF	RF	SST (Serum Separator)		
RPR	RPR	RPR	Gold (serum Separator)		
Rubella	RUB	RUB	Gold (serum Separator)		
Valproic Acid	VPA	VPA	Light Green (Lithium Heparin)		
Vancomycin	VANCO	VANCO	Light Green (Lithium Heparin)		
Urine Drugs of Abuse Screen ED and Inpatient	DAU	DRUGS	Urine		OPD order MISCE test
Blood Bank					
Type and Screen	TYHD	TYHD	Red Glass (No additive)		
Cultures					
Blood	BLC	BLC	Blood Culture bottle		
Chylamydia & GC DNA	CTNGQ	CTNGQ	Culture/urine		Culture or Urine
Throat Culture	THRC	THRC	culture		
MRSA-Nares	MRSAC	MRSAC	nasal culture		
Urine Culture	URNC	URNC	urine		
Vaginal R/O Group B Strep	GBSC	GBSC	vaginal culture		
Vaginal R/O Yeast + BV	BVSCRN	BVSCRN	vaginal culture		
Cdiff AG/TX with PCR reflex	CDIF	CDIFF2	stool		
Cryptosporidium AG	CRYPT	CRYPT	sterile stool specimen		
Stool Culture	STLC	STLC	sterile stool specimen		
Cyclospora	SPORA	SPORA	black stool container		Total fix
Giardia AG	GIARD	GIARD	plain stool specimen		
Ova & Parasites	OP	OP	sterile stool specimen		Total fix
Rapid Flu Screen	FLUAB	FLU	nasopharangeal swab		

TEST	McKesson	Lab	Container	Amount	Special Consideration
Hepatitis					
Hep A AB Total	HAVABQ	HAVABQ	SST (Serum Separator)		
Hep A IGM	HAM	HAM	SST (Serum Separator)		
Hep B Core AB Total	HBCAQ	HBCAQ	SST (Serum Separator)		
Hep B Surface AB	HBSB	HBSB	SST (Serum Separator)		
Hep B Surface AG	HBSG	HBSG	SST (Serum Separator)		
Hep C AB	HCB	HCAB	SST (Serum Separator)		
Hep C Virus RNA Quant	HCVRN	HCVRN	SST (Serum Separator)		
JR 2/5/2015					